


Dynalon Labware

Original Components
Proven Performance

See inside
to find out
more!



QWS4/120V 2.4 Liters/Hour Single Distillation Water Still

Introduction

Quickfit and Quartz launched their first ever Water Stills more than 60 years ago. Now Scilabware is proud to re-introduce the **Quickfit**® name into the water distillation market with the new **Quickfit**® QWS4 Water Still.

Manufactured in the UK, the QWS4 combines high performance with value for money. Other Water Stills may look similar, but only the QWS4 features a high-quality **Pyrex**® glass boiler and condenser.

Key Features

- Easy to install and easy to use
- Reliable operation – proven quality components
- High-quality distillate complies to ASTM and ISO standards
- Safe operation – heating element has two built-in safety thermostats
- Easy to clean – no glass dismantling required
- CE marked



Safety First

The chromium plated heating element is fitted with two thermostats (230V only) to protect the still from overheating in the event of water supply failure.

The thermostats can be reset to resume Water Still operation when the water supply has been restored.

Quality Distillate

The QWS4/120V will produce 2.4L/hr of pyrogen-free* distilled water from its reliable 3kW UK manufactured, chromium plated heating element.

Efficient Cooling


The highly efficient **Pyrex**® glass coil condenser, fitted with detachable plastic screw thread connectors, ensures that the distillate is cooled ready for use.



Easy Cleaning

No dismantling required. The **Pyrex**® glass boiler is fitted with a built-in funnel to allow the quick and easy addition of cleaning solutions.

Once the cleaning solution has done its job, simply open the GP **Rotaflo**® PTFE stopcock to drain the solution out of the boiler.



Technical Specification

Model. No.	540004-0120
Water Still Output	2.4 Liters / Hour
Distillate Quality:	
pH	5.0–6.5
Conductivity	3.0–4.0 μ S/cm
Resistivity	0.25–0.30megOhm-cm
Distillate Temperature	25–35°C
Water Supply Requirements:	
Water Flow Rate	1 litre/minute
Water Pressure	Minimum pressure 3 p.s.i. (0.207 bar)

*Care must **always** be taken to obtain pyrogen-free distilled water. The distillate **must** be tested prior to use.

Ordering Information

Product Code	Description
540004	Quickfit® 4Litres/Hour Single Distillation Water Still 230V 50/60Hz 3kW
540004-0120	Quickfit® 2.4 Litres/Hour Single Distillation Water Still 120V 50/60Hz 1.8kW
Spare Parts List:	
540014-0002	Pyrex® Boiler
540014-0003	Pyrex® Condenser
540014-0001	3kW Heating Element 230V 50/60Hz 3kW (with 2 Thermostats)
540014-0004	3kW Heating Element 120V 50/60Hz 1.8kW (with 1 Thermostat)
540014-0006	Hose Kit
Accessories:	
505634-2	Azlon® 10Litre Heavy Duty Polypropylene Aspirator Bottle
505634-3	Azlon® 20Litre Heavy Duty Polypropylene Aspirator Bottle

Pyrex® and Rotaflo® are registered trademarks. SciLabware registered user.

For more information visit www.dynalon.com



Dynalon® Labware

Email us with your questions loninfo@dyna-labware.com
Dynalab Corp. 175 Humboldt Street, Suite 300 Rochester, NY 14610 USA
Tel: 800.334.7585 | Fax: 585.334.0241