



**ViscoDispo**

The ViscoDispo disposable sampler can be used to take all layer and target point samples of high-viscosity media up to 100.000 mPas. The ergonomically shaped handle and the very stable drawbar make it easy to draw even viscous substances.

- Ergonomic handle
- No seals
- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-3320 | 500 mm | 100 mL | standard   | 20            |
| 853933-3321 | 500    | 100    | sterilized | 20            |
| 853933-3330 | 1000   | 200    | standard   | 20            |
| 853933-3331 | 1000   | 200    | sterilized | 20            |



**LiquiDispo**

The LiquiDispo disposable sampler can be used to take all layer and target point samples sample of liquids from highly liquid to slightly viscous. An ergonomic handle makes it easier to draw the sample: up to 200 ml can be drawn, as with a syringe.

- Ergonomic handle
- No seals
- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-1120 | 500 mm | 100 mL | standard   | 20            |
| 853933-1121 | 500    | 100    | sterilized | 20            |
| 853933-1130 | 1000   | 200    | standard   | 20            |
| 853933-1131 | 1000   | 200    | sterilized | 20            |

[www.dynalon.com](http://www.dynalon.com)



**DispoTube**

The ergonomic handle of the DispoTube makes it easier to use and facilitates a good grip even under difficult conditions, for example with gloves. Disposable sampler for taking samples of liquids.

- Ergonomic handle
- All-layer and target point samples possible
- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-6620 | 500 mm | 100 mL | standard   | 20            |
| 853933-6621 | 500    | 100    | sterilized | 20            |
| 853933-6630 | 1000   | 200    | standard   | 20            |
| 853933-6631 | 1000   | 200    | sterilized | 20            |



**DispoPipette**

Simple disposable sampler for taking samples of low viscosity media. If an opened DispoPipette is inserted into the liquid an all-layer sample is obtained. If it is closed while it is inserted up to the sampling point and then opened, a targetpoint sample is obtained. Very easy to operate. With two opening sizes for different viscosities.

- All-layer and target point samples possible
- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-5500 | 500 mm | 100 mL | standard   | 20            |
| 853933-5501 | 500    | 100    | sterilized | 20            |
| 853933-5531 | 1000   | 200    | standard   | 20            |
| 853933-5532 | 1000   | 200    | sterilized | 20            |

# LaboPlast® | SteriPlast® DispoSampler



- ▶ For disposable use
- ▶ Individually packaged
- ▶ Production, assembly and packaging in class 7 clean room (10,000)
- ▶ With EU foodstuffs and FDA approval
- ▶ Optionally sterilised by gamma rays





### DispoSampler – disposable sampler | for the highest standards in sampling

Disposable samplers – DispoSampler – are used to take samples of liquids, viscous media, powders and granulates. Target point, all-layer and cross-section samples can be taken, as well as small and large sample quantities. Unlike conventional reusable samplers, disposable samplers do not need to be cleaned after use.

- Cleaning reusable samplers is a demanding and often complicated task, with greasy, strongly adhesive or toxic substances. Sometimes the task is even impossible.
- Work processes do not have to be interrupted by cleaning
- Cross contamination is completely eliminated
- The cleaning method does not need to be validated in audits
- Time and money are saved

[www.dynalon.com](http://www.dynalon.com)



#### MultiDispo

Disposable zone samplers for taking samples from multiple layers of the sample material. Suitable for free-flowing powders and small-grain granulates. The representative all-layer sample can be emptied directly through the open handle into a sample container.

- Easy emptying through the opening in the handle
- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-4440 | 500 mm | 100 mL | standard   | 20            |
| 853933-4441 | 500    | 100    | sterilized | 20            |
| 853933-4450 | 1000   | 200    | standard   | 20            |
| 853933-4451 | 1000   | 200    | sterilized | 20            |



#### ZoneDispo

The ZoneDispo disposable all-layer sampler is used in sampling to collect a representative all-layer sample from three different layers. Suitable for free-flowing powders and small-grain granulates. The collected sample can be secured and transported in the closed sampler.

- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-4420 | 500 mm | 100 mL | standard   | 20            |
| 853933-4421 | 500    | 100    | sterilized | 20            |
| 853933-4430 | 1000   | 200    | standard   | 20            |
| 853933-4431 | 1000   | 200    | sterilized | 20            |



#### TargetDispo

TargetDispo is used to take a target point sample at a position in the sample material that has previously been precisely determined. The sampler pierces into the powder or granulate while closed. When it reaches the sample removal point it is opened and then closed again. In this way a precise target-point sample is obtained. The sample can simply be emptied through the open handle.

- Easy emptying through the opening in the handle
- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-4460 | 500 mm | 100 mL | standard   | 20            |
| 853933-4461 | 500    | 100    | sterilized | 20            |
| 853933-4470 | 1000   | 200    | standard   | 20            |
| 853933-4471 | 1000   | 200    | sterilized | 20            |

#### PowderDispo

PowderDispo is the ideal disposable sampler for cross-section samples of powdery or slightly cohesive substances. A complete all-layer sample can be taken from bags or smaller containers with the lance.

- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-5510 | 300 mm | 60 mL  | standard   | 20            |
| 853933-5511 | 300    | 60     | sterilized | 20            |
| 853933-5522 | 500    | 100    | standard   | 20            |
| 853933-5521 | 500    | 100    | sterilized | 20            |



#### DispoPicker

Disposable sampler for powders and small-grain bulk materials. Because the tip is closed, the sample material cannot trickle back out of the sample lance. DispoPicker is ideally suited for cross-section samples.

- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-5540 | 300 mm | 60 mL  | standard   | 20            |
| 853933-5541 | 300    | 60     | sterilized | 20            |
| 853933-5550 | 500    | 100    | standard   | 20            |
| 853933-5551 | 500    | 100    | sterilized | 20            |



#### DispoLance

The DispoLance disposable sampler can be used for taking large and small sample quantities quickly and easily, directly from bags. The sampler is suitable for free-flowing powders and small-grain granulates. The DispoLance sampler pierces into the container with the sturdy tip and the sample flows through the open inner tube directly into the container or bag that is held underneath.

- Variable sample quantity
- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Model      | Pack Quantity |
|-------------|--------|------------|---------------|
| 853933-5560 | 300 mm | standard   | 20            |
| 853933-5561 | 300    | sterilized | 20            |
| 853933-5570 | 500    | standard   | 20            |
| 853933-5571 | 500    | sterilized | 20            |



#### MicroDispo

Especially small sample quantities (10 ml) can be taken with MicroDispo. In pharmaceutical applications especially, often only a small quantity can or should be extracted for analysis. Using the disposable sampler, precise target point samples can be taken of free-flowing powders and small-grain bulk materials.

- Ergonomic handle
- HDPE
- Clean room manufactured
- Optionally sterilized by gamma rays

| Item        | Length | Volume | Model      | Pack Quantity |
|-------------|--------|--------|------------|---------------|
| 853933-2240 | 500 mm | 10 mL  | standard   | 20            |
| 853933-2241 | 500    | 10     | sterilized | 20            |
| 853933-2250 | 1000   | 10     | standard   | 20            |
| 853933-2251 | 1000   | 10     | sterilized | 20            |

