




# GHS Compliant Multi-lingual Wash Bottles



Dynalon  
Labware



- Multi-lingual wide mouth wash bottles featuring the patented DripLok® vapor venting valve that prevents solvent drips
- Indelibly printed in 4 languages (English, French, German and Spanish)
- Choose from four solvent labels, all including the new Globally Harmonized Symbols, NFPA diamond, risk phrases and safety information
- Fitted with color coded one-piece PP closures that maintain their shape in use
- Closures feature precision jet-tip aperture for fine control

Capacity mL	250	500	1000
Height mm	145	165	209
Diameter mm	60	75	90


	250 mL	500 mL
Printed	506485	506495
Acetone	-0001	-0001
Ethanol	-0002	-0002
Isopropanol	-0003	-0003
Methanol	-0004	-0004
Set/5*	-0005	-0005
Printed	506475	506475
Custom	-0001	-0002

506474-0003 1000mL Custom Wash Bottle

\*Sets include 1 of each printed bottle as well as the corresponding capacity custom write-on bottle.



# 500mL Capacity Wide Mouth Wash Bottles



- GHS Wide Mouth 500mL Wash Bottles feature the new Globally Harmonized Symbols, NFPA diamond, risk phrases and safety information
- Supplied in shrink wrapped pack of 5 wash bottles
- Closures with precision jet-tip aperture for fine control
- Fitted with white one-piece PP closures that maintain their shape in use

- A GHS Multilingual Wash Bottles, PP**
  - Non-venting and printed with Dichloromethane or Methyl ethyl ketone
  - Great for harsher chemicals and sterile environments
  - Indelibly printed in 4 languages (English, French, German and Spanish)
- B GHS Bleach Wash Bottle, Opaque LDPE**
  - Multi-lingual wide neck wash bottle
  - Opaque white LDPE for light sensitive sodium hypochlorite / bleach
  - Indelibly printed in three languages - English, French and Spanish
- C Autoclavable Wash Bottle, PP**
  - Autoclavable Wash Bottle at 121°C
  - Polypropylene 500mL bottle fitted with one-piece PP closure

Part No.	Printed	Capacity mL
506495-0015	Dichloromethane	500
506495-0016	Methyl ethyl ketone	500
506975	Sodium Hypochlorite	500
506985	Autoclavable	500

