


LTS Tips




GUIDELINES FOR USE

Light Touch Pipette Tips 85% Less Force And Higher Accuracy

LTS pipette tips users can achieve reliable and consistent fitting effect without excessive force. LTS pipette tips feature well-designed small sealing area and front stop point design. The front stop design allows you to know the exact position of the fitting.



Tip Ejection Forces



Please identify the LTS Pipette before purchasing

Multi-channel pipetting can be performed smoothly with the LTS tips system, ensuring consistent sample extraction across all channels. The front stop point design means that all channels are sealed perfectly at once, eliminating the need to shake or tighten the tips by hand.



**Dynalab®
Labware**




20ul Fit LTS Pipette, sterile

Part No.	With Filter	Low Retention	Package Type	Pack
RT020R	-	-	rack pack	96/rack, 50racks/case
RT020R-L	-	✓	rack pack	96/rack, 50racks/case
RT020RF	✓	-	rack pack	96/rack, 50racks/case
RT020RF-L	✓	✓	rack pack	96/rack, 50racks/case




200ul Fit LTS Pipette, sterile

Part No.	With Filter	Low Retention	Package Type	Pack
RT200R	-	-	rack pack	96/rack, 50racks/case
RT200R-L	-	✓	rack pack	96/rack, 50racks/case
RT200RF	✓	-	rack pack	96/rack, 50racks/case
RT200RF-L	✓	✓	rack pack	96/rack, 50racks/case



300ul Fit LTS Pipette, sterile

Part No.	With Filter	Low Retention	Package Type	Pack
RT300R	-	-	rack pack	96/rack, 50racks/case
RT300R-L	-	✓	rack pack	96/rack, 50racks/case
RT300RF	✓	-	rack pack	96/rack, 50racks/case
RT300RF-L	✓	✓	rack pack	96/rack, 50racks/case



1000ul Fit LTS Pipette, sterile

Part No.	With Filter	Low Retention	Package Type	Pack
RT1000R	-	-	rack pack	96/rack, 50racks/case
RT1000R-L	-	✓	rack pack	96/rack, 50racks/case
RT1000RF	✓	-	rack pack	96/rack, 50racks/case
RT1000RF-L	✓	✓	rack pack	96/rack, 50racks/case



1200ul Fit LTS Pipette, sterile

Part No.	With Filter	Low Retention	Package Type	Pack
RT1200R	-	-	rack pack	96/rack, 50racks/case
RT1200R-L	-	✓	rack pack	96/rack, 50racks/case
RT1200RF	✓	-	rack pack	96/rack, 50racks/case
RT1200RF-L	✓	✓	rack pack	96/rack, 50racks/case



Dynalon[®]
Labware

A CALIBRE SCIENTIFIC COMPANY

cs@dyna-labware.com
1 800-334-7585 | www.dynalon.com

