

Large Volume Tips



GUIDELINES FOR USE

100% pure virgin medical grade polypropylene.

Every lot of tips is tested and certified to be free of detectable Nuclease (Rnase/Dnase), Adenosine Tri-phosphate, and Endotoxin (pyrogen) contamination.

Large volume (5 ml, 10 ml) pipette tips can dispense large volume samples with the same accuracy as smaller volume tips.




Product Characteristics	Description
Sterilization grade	SAL 10-6.
Certified quality	Our products are made to ISO standards and we are audited lot to lot to ensure continuous quality.
Packing	ISO Certified.
Compatibility	Fit Eppendorf, Biohit, DLAB, Thermo.

5ml, clear, fit Eppendorf, Biohit, DLAB Pipette




Part No.	With Filter	Sterile	Pack
PT-5000R-E	<input type="checkbox"/>	<input checked="" type="checkbox"/>	24/rack, 30racks/case
PT-5000RF-E	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	24/rack, 30racks/case

10ml, clear, fit Eppendorf, Biohit, DLAB Pipette




Part No.	With Filter	Sterile	Pack
PT-10000R-E	<input type="checkbox"/>	<input checked="" type="checkbox"/>	24/rack, 40racks/case
PT-10000RF-E	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	24/rack, 40racks/case

5ml, clear, fit Thermo, Brand Pipette



Part No.	With Filter	Sterile	Pack
PT-5000R-T	<input type="checkbox"/>	<input checked="" type="checkbox"/>	40/rack, 40racks/case
PT-5000RF-T	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	40/rack, 40racks/case

10ml, clear, fit Thermo, Brand Pipette



Part No.	With Filter	Sterile	Pack
PT-10000R-T	<input type="checkbox"/>	<input checked="" type="checkbox"/>	24/rack, 40racks/case
PT-10000RF-T	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	24/rack, 40racks/case



**Dynamalab®
Labware**



Compatibility List

Product No.	Pipette Volume	Pipette Brand									
		Eppendorf		Thermo		Biohit		DLAB		Brand	
		Research plus 0.5-5ml	Research plus 1-10ml	Finnpipette F3 0.5-5ml	Finnpipette F3 1-10ml	Proline Plus 500-5000ul	Proline Plus 1-10ml	HiPette 1000-5000ul	HiPette 2-10ml	Transferpette® S 500-5000ul	Transferpette® S 1000-10000ul
PT-5000R-E	500-5000ul	✓				✓		✓			
PT-5000RF-E	500-5000ul	✓				✓		✓			
PT-10000R-E	1000-10000ul		✓				✓		✓		
PT-10000RF-E	1000-10000ul		✓				✓		✓		
PT-5000R-T	500-5000ul			✓						✓	
PT-5000RF-T	500-5000ul			✓						✓	
PT-10000R-T	1000-10000ul				✓						✓
PT-10000RF-T	1000-10000ul				✓						✓



Dynalabware

A CALIBRE SCIENTIFIC COMPANY

cs@dyna-labware.com
1 800-334-7585 | www.dynalab.com

