

Serological Pipettes



GUIDELINES FOR USE

Product Specification

Our rigorous quality specifications and convenient functionality make our serological pipettes an excellent choice for routine laboratory work and cell studies. Bold, clear bidirectional scale makes liquid capacity and sample addition easy to read. A variety of specifications and packaging options are available. Certified scale accuracy and precision calibration ensures $\pm 2\%$ accuracy in the 1mL and 2mL total volume range and $\pm 1\%$ accuracy in the 5mL, 10mL, 25mL, 50mL total volume range.

- All of our serological pipettes are tested by an independent laboratory and certified to be RNase\ DNase\PCR inhibitors\human DNA\ATP\Pyrogen free.
- ISO 13485:2016 Quality control system.
- Materials meet the requirements of USP Class VI.
- Irradiation sterilization, sterility grade SAL 10-6.
- Universal design makes pipette suitable for most brands of assistance devices.
- Registered abroad as in vitro diagnostic supplies.

Material Information

High transparency polystyrene meets USP Class VI medical grade material requirements.

Product Characteristics

- Filter plugs prevent aerosol or liquid contamination of the pipette unit and further reduce the risk of contamination in the sample room.
- Color coded serum pipette and packaging for easy selection of volume.
- Serological pipette has high sterility guarantee of paper/plastic or plastic easy to tear independent packaging.
- Bulk package option available reducing packaging waste.
- Validity period: 3 years from the production date.



Dynalon[®]
Labware

A CALIBRE SCIENTIFIC COMPANY



● 1ml Serological Pipette, PS, Sterile


Part No.	Package Type	Pack
SPTS01	Individual pack	50/bag, 20bags/case
SPTB01	bulk pack	50/bag, 20bags/case

● 5ml Serological Pipette, PS, Sterile

Part No.	Package Type	Pack
SPTS05	Individual pack	50/bag, 10bags/case
SPTB05	bulk pack	50/bag, 10bags/case

● 10ml Serological Pipette, PS, Sterile

Part No.	Pcs	Package Type	Pack
SPTS10	one pcs	Individual pack	50/bag, 10bags/case
SPTB10	one pcs	bulk pack	50/bag, 10bags/case
SPTS10A	two pcs	Individual pack	50/bag, 10bags/case
SPTB10A	two pcs	bulk pack	50/bag, 10bags/case



● 2ml Serological Pipette, PS, Sterile


Part No.	Package Type	Pack
SPTS02	Individual pack	50/bag, 20bags/case
SPTB02	bulk pack	50/bag, 20bags/case

● 25ml Serological Pipette, PS, Sterile

Part No.	Package Type	Pack
SPTS25	Individual pack	25/bag, 10bags/case
SPTB25	bulk pack	25/bag, 10bags/case

● 50ml Serological Pipette, PS, Sterile

Part No.	Package Type	Pack
SPTS50	Individual pack	25/bag, 8bags/case
SPTB50	bulk pack	25/bag, 8bags/case



Size (ml)	Appearance Size												
	A	B	C	D	L1	ID1	OD1	L3	ID3	OD3	L2	ID2	OD2
1	266±2	2.3±0.3	4.6±0.2	4.6±0.2	-	1.6±0.3	2.3±0.3	-	2.8±0.2	4.6±0.2	-	2.8±0.2	4.6±0.2
2	270±2	2.7±0.3	6.3±0.2	6.3±0.2	-	1.7±0.3	2.7±0.3	-	4.1±0.2	6.3±0.2	-	4.1±0.2	6.3±0.2
5	342±2	3.2±0.3	8.1±0.2	8.1±0.2	-	2.2±0.3	3.2±0.3	-	6±0.2	8.1±0.2	-	6±0.2	8.1±0.2
10	344±2	3.6±0.15	7.9±0.2	10±0.2	37	1.5±0.1	3.6±0.15	19.2±0.2	5.5±0.2	7.9±0.2	285	8.2±0.2	10±0.2
25	313±2	5±0.2	8.1±0.2	15±0.2	35	2.5±0.2	5±0.2	21.7±0.2	5.2±0.2	8.1±0.2	255	13±0.2	15±0.2
50	353±2	5.6±0.2	8.1±0.2	18±0.2	35	3±0.2	5.6±0.2	21.7±0.2	5.2±0.2	8.1±0.2	295	16±0.2	18±0.2



Dynalon[®]
Labware

A CALIBRE SCIENTIFIC COMPANY

cs@dyna-labware.com
1 800-334-7585 | www.dynalon.com

