

## Technical Datasheet

|                                 |  |              |
|---------------------------------|--|--------------|
| <b>Article - No:</b>            | TDS-23-057   | Revision: A0 |
| <b>Remarks</b>                  | single use only  |              |
| <b>Specification</b>            |  |              |
| <b>Product name</b>             | CL-D060, 60mm, Cell Culture Dish                       |              |
| <b>Production conditions</b>    | Class 100,00 clean room production                     |              |
|                                 | No release agent is used                               |              |
| <b>Classification</b>           | GMP Compliant  |              |
| <b>Dimensions</b>               |  |              |
| <b>Colour</b>                   | Clear  |              |
| <b>Nominal Capacity</b>         | 60mm   |              |
| <b>Material Characteristics</b> | polystyrene (PS), premium surface                      |              |
| <b>Sterile</b>                  | Yes  |              |
| <b>Product Characteristics</b>  | biologically inert                                     |              |
|                                 | physiologically harmless                               |              |
|                                 | no heavy metals according to EC/94/62                  |              |
| <b>Purity Grade</b>             | free of RNase and DNase                                |              |
|                                 | free of pyrogens                                       |              |
|                                 | free of PCR inhibitors                                 |              |
|                                 | free of human DNA                                      |              |
|                                 | free of endotoxins                                     |              |
|                                 | free of ATP  |              |
| <b>Shelf - Life</b>             | 36 months at controlled storage conditions             |              |
| <b>Packaging</b>                | pcs per bag: 10  |              |
|                                 | bag per case: 50                                       |              |
| <b>Storage</b>                  | Storage Temperature: 20 ° C-25 ° C, humidity 40% - 60% |              |
|                                 | Keep clean, ventilated and out of direct sunlight      |              |