


Technical Datasheet

Article - No:	TDS-23-056	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-DT035, 35mm, Cell Culture Dish, TC Treated	
Production conditions	Class 100,00 clean room production	
	No release agent is used	
Classification	GMP Compliant	
Dimensions	 <p>The drawing shows a side view of the dish with a diameter of 38.80 mm. The height of the dish is 12.45 mm, and the height of the rim is 6.50 mm. The top view shows a circular dish with a diameter of 38.80 mm and a central hole. The drawing is labeled with 'A' at the top and bottom, and 'Ø38.80' for the diameter.</p>	
Colour	Clear	
Nominal Capacity	35mm	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 20	
	bag per case: 25	
Storage	Storage Temperature: 20 ° C-25 ° C, humidity 40% - 60%	
	Keep clean, ventilated and out of direct sunlight	