

Technical Datasheet

Article - No:	TDS-23-060	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-DT100, 100mm, Cell Culture Dish, TC Treated	
Production conditions	Class 100,00 clean room production	
	No release agent is used	
Classification	GMP Compliant	
Dimensions		
Colour	Clear	
Nominal Capacity	100mm	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 10	
	bag per case: 30	
Storage	Storage Temperature: 20 ° C-25 ° C, humidity 40% - 60%	
	Keep clean, ventilated and out of direct sunlight	