


Technical Datasheet

Article - No:	TDS-23-003	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-F025V, 25cm ² , Vent Cap, Cell Culture Flask	
Production conditions	Class 100,00 clean room production	
	No release agent is used	
Classification	GMP Compliant	
Dimensions	 <p>The drawing shows a side view of a cell culture flask. It has a rectangular body with rounded corners and a tapered neck leading to a cylindrical vent cap. Three dimensions are indicated with arrows: a total length of 97.45 units, a length of 76.24 units from the left edge to the start of the vent cap, and a height of 49.44 units.</p>	
Colour	Clear	
Graduation	printed rings	
Nominal Capacity	25cm ²	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 10	
	bag per case: 20	
Storage	Storage Temperature: 20 ° C-25 ° C, humidity 40% -60%	
	Keep clean, ventilated and out of direct sunlight	