


## Technical Datasheet

<b>Article - No:</b>	TDS-23-005	Revision: A0
<b>Remarks</b>	single use only	
<b>Specification</b>		
<b>Product name</b>	CL-F075P, 75cm <sup>2</sup> , Plug Cap, Cell Culture Flask	
<b>Production conditions</b>	Class 100,00 clean room production	
	No release agent is used	
<b>Classification</b>	GMP Compliant	
<b>Dimensions</b>	 <p>The drawing shows a side view of a cell culture flask. It has a rectangular body with rounded corners, a tapered neck, and a cylindrical plug cap. Three dimensions are indicated with arrows: a total length of 154.7 units, a length of 129.2 units from the left edge to the start of the neck, and a height of 81.1 units.</p>	
<b>Colour</b>	Clear	
<b>Graduation</b>	printed rings	
<b>Nominal Capacity</b>	75cm <sup>2</sup>	
<b>Material Characteristics</b>	polystyrene (PS), premium surface	
<b>Sterile</b>	Yes	
<b>Product Characteristics</b>	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
<b>Shelf - Life</b>	36 months at controlled storage conditions	
<b>Packaging</b>	pcs per bag: 5	
	bag per case: 20	
<b>Storage</b>	Storage Temperature: 20 ° C-25 ° C, humidity 40% -60%	
	Keep clean, ventilated and out of direct sunlight	