


Technical Datasheet

| | | |
|---------------------------------|---|--------------|
| Article - No: | TDS-23-002 | Revision: A0 |
| Remarks | single use only | |
| Specification | | |
| Product name | CL-FT025P, 25cm ² , Plug Cap, Cell Culture Flask, TC Treated | |
| Production conditions | Class 100,00 clean room production | |
| | No release agent is used | |
| Classification | GMP Compliant | |
| Dimensions |  <p>The drawing shows a side view of a cell culture flask. It has a rectangular body with rounded corners, a tapered neck, and a cylindrical cap. Dimension lines indicate: a total length of 97.45 units, a length of 76.24 units from the left edge to the start of the neck, and a height of 49.44 units.</p> | |
| Colour | Clear | |
| Graduation | printed rings | |
| Nominal Capacity | 25cm ² | |
| Material Characteristics | polystyrene (PS), premium surface | |
| Sterile | Yes | |
| Product Characteristics | biologically inert | |
| | physiologically harmless | |
| | no heavy metals according to EC/94/62 | |
| Purity Grade | free of RNase and DNase | |
| | free of pyrogens | |
| | free of PCR inhibitors | |
| | free of human DNA | |
| | free of endotoxins | |
| | free of ATP | |
| Shelf - Life | 36 months at controlled storage conditions | |
| Packaging | pcs per bag: 10 | |
| | bag per case: 20 | |
| Storage | Storage Temperature: 20 ° C-25 ° C, humidity 40% - 60% | |
| | Keep clean, ventilated and out of direct sunlight | |