


## Technical Datasheet

<b>Article - No:</b>	TDS-23-008	Revision: A0
<b>Remarks</b>	single use only	
<b>Specification</b>		
<b>Product name</b>	CL-FT075V, 75cm <sup>2</sup> , Vent Cap, Cell Culture Flask, TC Treated	
<b>Production conditions</b>	Class 100,00 clean room production	
	No release agent is used	
<b>Classification</b>	GMP Compliant	
<b>Dimensions</b>	 <p>The drawing shows a side view of a cell culture flask with a vent cap. The total length of the flask is 154.7 units. The length from the start of the neck to the center of the vent cap is 129.2 units. The height of the flask at its widest point is 81.1 units.</p>	
<b>Colour</b>	Clear	
<b>Graduation</b>	printed rings	
<b>Nominal Capacity</b>	75cm <sup>2</sup>	
<b>Material Characteristics</b>	polystyrene (PS), premium surface	
<b>Sterile</b>	Yes	
<b>Product Characteristics</b>	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
<b>Shelf - Life</b>	36 months at controlled storage conditions	
<b>Packaging</b>	pcs per bag: 5	
	bag per case: 20	
<b>Storage</b>	Storage Temperature: 20 ° C-25 ° C, humidity 40% -60%	
	Keep clean, ventilated and out of direct sunlight	