


Technical Datasheet

Article - No:	TDS-23-0010	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-FT175P, 175cm ² , Plug Cap, Cell Culture Flask, TC Treated	
Production conditions	Class 100,00 clean room production	
	No release agent is used	
Classification	GMP Compliant	
Dimensions	 <p>The drawing shows a side view of a cell culture flask with a plug cap. The main body length is 178.4 units, the total length including the neck and cap is 228.2 units, and the height of the main body is 115.6 units.</p>	
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	175cm ²	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 5	
	bag per case: 10	
Storage	Storage Temperature: 20 ° C-25 ° C, humidity 40% -60%	
	Keep clean, ventilated and out of direct sunlight	