

Technical Datasheet

Article - No:	PG-QA-R025-46	Revision: A0
Remarks	single use only	
Specification		
Product name	CT15R, 15ml Centrifuge tube	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVDD	
Dimensions		
Colour	Clear tube with Orange cap	
Graduation	Printing	
Nominal Capacity	15ml	
Material Characteristics	polypropylene (PP), premium surface	
Performance parameters	Storage Temperature: -80°C~121°C	
	Centrifugal force: 15000xg	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 25	
	rack per case: 20	
Storage	dry condition: 40%-60% rel. humidity	