


## Technical Datasheet

<b>Article - No:</b>	PG-QA-R025-09-004	Revision: A0																		
<b>Remarks</b>	single use only																			
<b>Specification</b>																				
<b>Product name</b>	DWPC96, KingFisher Plate Combo																			
<b>Production conditions</b>	Class 100,000 clean room production																			
	No release agent is used																			
	Motor molding																			
<b>Classification</b>	CE/IVD																			
<b>Dimensions</b>	 <table border="1" data-bbox="1053 728 1524 817"> <thead> <tr> <th>Part No.</th> <th>Description</th> <th>A ±0.3</th> <th>B ±0.2</th> <th>C ±0.2</th> <th>D ±0.2</th> <th>E ±0.2</th> <th>F ±0.1</th> <th>G ±0.2</th> </tr> </thead> <tbody> <tr> <td>DWPC96</td> <td>KingFisher Plate Combo</td> <td>127.50</td> <td>99.00</td> <td>85.30</td> <td>63.00</td> <td>109.70</td> <td>5.00</td> <td>44.90</td> </tr> </tbody> </table>		Part No.	Description	A ±0.3	B ±0.2	C ±0.2	D ±0.2	E ±0.2	F ±0.1	G ±0.2	DWPC96	KingFisher Plate Combo	127.50	99.00	85.30	63.00	109.70	5.00	44.90
	Part No.	Description	A ±0.3	B ±0.2	C ±0.2	D ±0.2	E ±0.2	F ±0.1	G ±0.2											
DWPC96	KingFisher Plate Combo	127.50	99.00	85.30	63.00	109.70	5.00	44.90												
<b>Colour</b>	Clear																			
<b>Bottle Shape</b>	V type																			
<b>Material Characteristics</b>	polypropylene (PP), premium surface																			
<b>Product Characteristics</b>	biologically inert																			
	physiologically harmless																			
	no heavy metals according to EC/94/62																			
<b>Purity Grade</b>	free of RNase and DNase																			
	free of endotoxins																			
<b>Shelf - Life</b>	36 months at controlled storage conditions																			
<b>Packaging</b>	pcs per bag: 10																			
	bag per case: 10																			
<b>Storage</b>	Storage Temperature: 2 ° C-40 ° C, humidity 40% -80%																			
	Keep clean, ventilated and out of direct sunlight																			