


Technical Datasheet

Article - No:	PG-QA-R025-09-002	Revision: A0																				
Remarks	single use only																					
Specification																						
Product name	DWPS20U, 2.2ml deep well plate, Square well U bottom																					
Production conditions	Class 100,000 clean room production																					
	No release agent is used																					
	Motor molding																					
Classification	CE/IVD																					
Dimensions	 <table border="1" data-bbox="895 808 1549 907"> <thead> <tr> <th>Part No.</th> <th>Description</th> <th>A ±0.2</th> <th>B ±0.2</th> <th>C ±0.2</th> <th>D ±0.1</th> <th>E ±0.3</th> <th>F ±0.05</th> <th>G ±0.2</th> <th>H +0.4/-0.2</th> </tr> </thead> <tbody> <tr> <td>DWPS20U</td> <td>2.2ml square well plate, round well U bottom</td> <td>127.30</td> <td>85.30</td> <td>41.60</td> <td>0.60</td> <td>2.95</td> <td>8.20</td> <td>123.13</td> <td>80.93</td> </tr> </tbody> </table>		Part No.	Description	A ±0.2	B ±0.2	C ±0.2	D ±0.1	E ±0.3	F ±0.05	G ±0.2	H +0.4/-0.2	DWPS20U	2.2ml square well plate, round well U bottom	127.30	85.30	41.60	0.60	2.95	8.20	123.13	80.93
Part No.	Description	A ±0.2	B ±0.2	C ±0.2	D ±0.1	E ±0.3	F ±0.05	G ±0.2	H +0.4/-0.2													
DWPS20U	2.2ml square well plate, round well U bottom	127.30	85.30	41.60	0.60	2.95	8.20	123.13	80.93													
Colour	Clear																					
Bottle Shape	U type																					
Nominal Capacity	2.2ml																					
Material Characteristics	polypropylene (PP), premium surface																					
Product Characteristics	biologically inert																					
	physiologically harmless																					
	no heavy metals according to EC/94/62																					
Purity Grade	free of RNase and DNase																					
	free of endotoxins																					
Shelf - Life	36 months at controlled storage conditions																					
Packaging	pcs per bag: 10																					
	bag per case: 10																					
Storage	Storage Temperature: 2 ° C-40 ° C, humidity 40% -80%																					
	Keep clean, ventilated and out of direct sunlight																					