

Technical Datasheet

| | | |
|---|---|--------------|
| Article - No: | TDS-PD150H | Revision: A0 |
| Remarks | single use only | |
| Specification | | |
| Product name | PD150H, 150mm Petri dish, 20mm height | |
| Production conditions | Class 100,000 clean room production | |
| | No release agent is used | |
| Classification | CE/IVD | |
| Dimensions | <p>The drawing shows four views of the petri dish components: a top view of the cap with a diameter of $\phi 143.6$, a side view of the cap with a height of 20.2 ± 0.3, a top view of the bottom with an inner diameter of 138.2 ± 0.5, and a side view of the bottom with a height of 20.2 ± 0.3.</p> | |
| Colour | Clear | |
| Graduation | N/A | |
| Nominal Capacity | 150mm diameter | |
| Material Characteristics | polystyrene (PS), premium surface | |
| Sterile | EO sterilization | |
| Product Characteristics | biologically inert | |
| | no heavy metals according to EC/94/62 | |
| Purity Grade | free of RNase and DNase | |
| | free of pyrogens | |
| | free of PCR inhibitors | |
| | free of human DNA | |
| | free of endotoxins | |
| Shelf - Life | free of ATP | |
| | 36 months at controlled storage conditions | |
| Packaging | pcs per bag: 20 | |
| | bag per case: 6 | |
| Storage | Storage Temperature: 20 ° C-25 ° C, humidity 40% -60% | |
| | Keep clean, ventilated and out of direct sunlight | |
| Approved by PakGent quality department | | |