

# Technical Datasheet

<b>Article - No:</b>	TDS-23-063	Revision: A0																							
<b>Remarks</b>	single use only																								
<b>Specification</b>																									
<b>Product name</b>	PD35, 35mm, Petri Dish																								
<b>Production conditions</b>	Class 100,00 clean room production																								
	No release agent is used																								
<b>Classification</b>	GMP Compliant																								
<b>Dimensions</b>																									
	<table border="1" style="display: inline-table; margin-right: 20px;"> <thead> <tr> <th>A</th> <th>B</th> <th>C</th> <th>D</th> <th>E</th> <th>F</th> </tr> </thead> <tbody> <tr> <td><math>\varphi 36.5 \pm 0.02</math></td> <td><math>\varphi 32.1</math></td> <td><math>0.8 \pm 0.02</math></td> <td><math>0.3 \pm 0.02</math></td> <td><math>11.0 \pm 0.02</math></td> <td><math>9.8 \pm 0.02</math></td> </tr> </tbody> </table> <table border="1" style="display: inline-table;"> <thead> <tr> <th>A</th> <th>B</th> <th>C</th> <th>D</th> <th>E</th> <th>F</th> </tr> </thead> <tbody> <tr> <td><math>\varphi 39.4 \pm 0.02</math></td> <td><math>\varphi 33.5 \pm 0.02</math></td> <td><math>8.2 \pm 0.02</math></td> <td>7.60</td> <td>4.00°</td> <td>R1.00</td> </tr> </tbody> </table>	A	B	C	D	E	F	$\varphi 36.5 \pm 0.02$	$\varphi 32.1$	$0.8 \pm 0.02$	$0.3 \pm 0.02$	$11.0 \pm 0.02$	$9.8 \pm 0.02$	A	B	C	D	E	F	$\varphi 39.4 \pm 0.02$	$\varphi 33.5 \pm 0.02$	$8.2 \pm 0.02$	7.60	4.00°	R1.00
A	B	C	D	E	F																				
$\varphi 36.5 \pm 0.02$	$\varphi 32.1$	$0.8 \pm 0.02$	$0.3 \pm 0.02$	$11.0 \pm 0.02$	$9.8 \pm 0.02$																				
A	B	C	D	E	F																				
$\varphi 39.4 \pm 0.02$	$\varphi 33.5 \pm 0.02$	$8.2 \pm 0.02$	7.60	4.00°	R1.00																				
<b>Colour</b>	Clear																								
<b>Nominal Capacity</b>	35mm																								
<b>Material Characteristics</b>	polystyrene (PS), premium surface																								
<b>Sterile</b>	Yes																								
<b>Product Characteristics</b>	biologically inert																								
	physiologically harmless																								
	no heavy metals according to EC/94/62																								
<b>Purity Grade</b>	free of RNase and DNase																								
	free of pyrogens																								
	free of PCR inhibitors																								
	free of human DNA																								
	free of endotoxins																								
	free of ATP																								
<b>Shelf - Life</b>	36 months at controlled storage conditions																								
<b>Packaging</b>	pcs per bag: 10																								
	bag per case: 230																								
<b>Storage</b>	Storage Temperature: 20 ° C-25 ° C, humidity 40% -60%																								
	Keep clean, ventilated and out of direct sunlight																								