

Technical Datasheet

Article - No:	TDS-23-064	Revision: A0																							
Remarks	single use only																								
Specification																									
Product name	PD60, 60mm, Petri Dish																								
Production conditions	Class 100,00 clean room production																								
	No release agent is used																								
Classification	GMP Compliant																								
Dimensions																									
	<table border="1"> <thead> <tr> <th>A</th> <th>B</th> <th>C</th> <th>D</th> <th>E</th> <th>F</th> <th>A</th> <th>B</th> <th>C</th> <th>D</th> <th>E</th> <th>F</th> </tr> </thead> <tbody> <tr> <td>φ36.5±0.02</td> <td>φ32.1</td> <td>0.8±0.02</td> <td>0.3±0.02</td> <td>11.0±0.02</td> <td>9.8±0.02</td> <td>φ39.4±0.02</td> <td>φ33.5±0.02</td> <td>8.2±0.02</td> <td>7.60</td> <td>4.00°</td> <td>RL00</td> </tr> </tbody> </table>		A	B	C	D	E	F	A	B	C	D	E	F	φ36.5±0.02	φ32.1	0.8±0.02	0.3±0.02	11.0±0.02	9.8±0.02	φ39.4±0.02	φ33.5±0.02	8.2±0.02	7.60	4.00°
A	B	C	D	E	F	A	B	C	D	E	F														
φ36.5±0.02	φ32.1	0.8±0.02	0.3±0.02	11.0±0.02	9.8±0.02	φ39.4±0.02	φ33.5±0.02	8.2±0.02	7.60	4.00°	RL00														
Colour	Clear																								
Nominal Capacity	60mm																								
Material Characteristics	polystyrene (PS), premium surface																								
Sterile	Yes																								
Product Characteristics	biologically inert																								
	physiologically harmless																								
	no heavy metals according to EC/94/62																								
Purity Grade	free of RNase and DNase																								
	free of pyrogens																								
	free of PCR inhibitors																								
	free of human DNA																								
	free of endotoxins																								
	free of ATP																								
Shelf - Life	36 months at controlled storage conditions																								
Packaging	pcs per bag: 20																								
	bag per case: 50																								
Storage	Storage Temperature: 20 ° C-25 ° C, humidity 40% - 60%																								
	Keep clean, ventilated and out of direct sunlight																								