


Technical Datasheet

Article - No:	PG-QA-R025-01-0075	Revision: A0
Remarks	single use only	
Specification		
Product name	PT-10000R-E, 10ML Pipette tips	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	clear	
Graduation	moulded rings	
Nominal Capacity	10ML	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 24	
	rack per case: 40	
Storage	Storage Temperature: 2 ° C-40 ° C, humidity 40% -80%	
	Keep clean, ventilated and out of direct sunlight	