


Technical Datasheet

Article - No:	TDS-SMB1000G	Revision: A0
Remarks	single use only	
Specification		
Product name	SMB1000G, 1000ml serum bottle, PETG	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVDD	
Dimensions		
Colour	Clear	
Nominal Capacity	1000ml	
Material Characteristics	bottle:PETG, surface	
	cap: High Density Polyethylene(HDPE)	
Performance parameters	Negative pressure: -70kpa~-90kpa	
	Ethanol loss rate: ≤0.5%	
	water loss rate: ≤0.2%	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 12	
	bag per case: 2	
Storage	dry condition: 40%-60% rel. humidity	