


## Technical Datasheet

<b>Article - No:</b>	TDS-SMB125G	Revision: A0
<b>Remarks</b>	single use only	
Specification		
<b>Product name</b>	SMB125G, 125ml serum bottle, PETG	
<b>Production conditions</b>	Class 100,000 clean room production	
	No release agent is used	
<b>Classification</b>	CE/IVDD	
<b>Dimensions</b>	 <p style="text-align: center;"> <math>\varnothing 28.80 \pm 0.1</math>            2.10            4.20         </p>	
<b>Colour</b>	Clear	
<b>Nominal Capacity</b>	125ml	
<b>Material Characteristics</b>	bottle: PETG, surface	
	cap: High Density Polyethylene (HDPE)	
<b>Performance parameters</b>	Negative pressure: -70kpa~-90kpa	
	Ethanol loss rate: $\leq 0.5\%$	
	water loss rate: $\leq 0.2\%$	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
<b>Shelf - Life</b>	36 months at controlled storage conditions	
<b>Packaging</b>	pcs per bga: 24	
	bag per case: 4	
<b>Storage</b>	dry condition: 40%-60% rel. humidity	