

Technical Datasheet

Article - No:	PG-QA-R025-01-0092	Revision: A0
Remarks	single use only	
Specification		
Product name	UFPT-REF-200-L, 200uL Pipette tips, Refull pack, Low retention	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p>Technical drawing of a pipette tip showing dimensions A through L. Dimension A is the tip length, B is the total length, C is the length to the start of the graduation, D is the length of the graduation, E is the length of the body, F is the diameter of the body, G is the diameter of the tip, H is the diameter of the tip, I is the diameter of the tip, and L is the diameter of the tip.</p>	
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	1-200UL	
Sterile	Yes	
Material Characteristics	polypropylene(PP),premium surface	
	pack:Polyethylene terephthalate(PET)	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per pack: 12	
	pack per case: 10	
	pcs per case: 11520	
Storage	Storage Temperature: 2 ° C-40 ° C, humidity 40% -80%	
	Keep clean, ventilated and out of direct sunlight	