

Technical Datasheet

Article - No:	PG-QA-R025-01-0091	Revision: A0
Remarks	single use only	
Specification		
Product name	UFPT-REF-200, 200uL Pipette tips, Refull pack	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p>Technical drawing of a pipette tip with dimensions:</p> <ul style="list-style-type: none"> A: Tip length B: Tip diameter at the top C: Total length: 50.6 ± 0.2 D: Length to the start of the rings: 28.00 ± 0.4 E: Length of the rings: 16.10 F: Diameter of the rings: $\Phi 6.30$ G: Diameter of the tip at the bottom: $\Phi 7.4 \pm 0.05$ H: Diameter of the narrow part: $\Phi 0.6$ I: Diameter of the wide part: $\Phi 0.8$ L: Diameter of the tip at the bottom 	
Colour	Yellow	
Graduation	moulded rings	
Nominal Capacity	1-200UL	
Sterile	Yes	
Material Characteristics	polypropylene(PP),premium surface	
	pack:Polyethylene terephthalate(PET)	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per pack: 12	
	pack per case: 10	
	pcs per case: 11520	
Storage	Storage Temperature: 2 ° C-40 ° C, humidity 40% -80%	
	Keep clean, ventilated and out of direct sunlight	