


Technical Datasheet

Article - No:	TDS-WMPB1000	Revision: A0																							
Remarks	single use only																								
Specification																									
Product name	WMPB1000,1000ml Packaging Bottle,HDPE																								
Production conditions	Class 100,000 clean room production																								
	No release agent is used																								
Classification	CE/IVDD																								
Dimensions	<div></div> <table><tr><th>location name</th><th>material</th><th>O. D. D1 (mm)</th><th>Height H1 (mm)</th><th>I. D. D2 (mm)</th><th>Cap O.D D3 (mm)</th><th>Cap H H2 (mm)</th><th>Total H (mm)</th></tr><tr><td rowspan="2">1000ml</td><td>PP</td><td>φ91.9±0.5</td><td>193.9±0.5</td><td>φ43.7±0.4</td><td>φ57.5±0.4</td><td>16.8±0.4</td><td>196.0±1.0</td></tr><tr><td>PE</td><td>φ92.2±0.5</td><td>193.9±0.5</td><td>φ43.7±0.4</td><td>φ57.5±0.4</td><td>16.8±0.4</td><td>196.0±1.0</td></tr></table>		location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)	1000ml	PP	φ91.9±0.5	193.9±0.5	φ43.7±0.4	φ57.5±0.4	16.8±0.4	196.0±1.0	PE	φ92.2±0.5	193.9±0.5	φ43.7±0.4	φ57.5±0.4	16.8±0.4	196.0±1.0
location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	Total H (mm)																		
1000ml	PP	φ91.9±0.5	193.9±0.5	φ43.7±0.4	φ57.5±0.4	16.8±0.4	196.0±1.0																		
	PE	φ92.2±0.5	193.9±0.5	φ43.7±0.4	φ57.5±0.4	16.8±0.4	196.0±1.0																		
Colour	Clear																								
Nominal Capacity	1000ml																								
Material Characteristics	bottle:High Density Polyethylene(HDPE)																								
	cap:polypropylene (PP)																								
Performance parameters	Negative pressure:-70kpa~-90kpa																								
	Ethanol loss rate:≤0.5%																								
	water loss rate:≤0.2%																								
Purity Grade	free of RNase and DNase																								
	free of pyrogen																								
	free of PCR inhibitors																								
	free of human DNA																								
	free of endotoxins																								
Shelf - Life	36 months at controlled storage conditions																								
Packaging	pcs per case:50																								
Storage	dry condition:40%-60% rel. humidity																								
Approved by quality department																									