




## Ideal for Danio rerio breeding Dura-Cross™ Zebrafish Tank



Dynalon  
Labware




Dura-Cross Zebrafish breeding tank is manufactured from transparent and autoclavable food grade polycarbonate (PC), it is fully impact resistant and will not chip or break. The tank consists of a tapered tank body, a slotted insert and lid, and is stackable to save space.

175 Humboldt Street, Suite 300  
Rochester, NY 14610  
Telephone 800-334-7585  
Fax 585-334-0241  
[www.dynalon.com](http://www.dynalon.com)




# Dura-Cross™ Zebrafish Breeding Tank Optional Dividers



Dura-Cross easily and effectively optimizes the breeding of zebra fish.

Features and accessories include:

- Manufactured from food grade polycarbonate - virtually indestructible thereby minimizing replacement costs. Glass clear for maximum viewing of fish and eggs.
- Unique slotted insert ensures eggs are separated from adults, maximizing yield; the removable insert allows adult fish to be quickly and safely returned to the main habitat tank.
- Molded lid keeps the tank and fish environment clean and is designed for easy stacking. Two 1/4" holes are able to take air tubing.
- Optional Dividers made of PC are available in diagonal and straight formats to keep males and females separate as required.
- Optional Carrying Tray made of PC allows for safe transport of tanks; tray holds 5 tanks in a single layer or 10 tanks stacked.
- Tanks and accessories are dishwasher safe and autoclavable, allowing for easy cleaning.
- Tank dimensions: 8 x 4 x 3" (20.3 x 10.2 x 7.6cm)



Catalog Number	Description	Price
143354-0000	Zebra Fish Tank Complete	\$40.85
143354-0002	Zebra Tank Divider - Horizontal	\$5.28
143354-0004	Zebra Tank Divider - Diagonal	\$5.55
143354-0006	Zebra Tank Carrier*	\$52.84

\*Holds 5 to 10 Zebra Fish Tanks

**Minimum order quantities required.**

Contact customer service [CS@dyna-labware.com](mailto:CS@dyna-labware.com) for further details.

