


Technical Datasheet

Article - No:	TDS-23-016	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-PT012, 12-well, Cell Culture Plate, TC Treated	
Production conditions	Class 100,00 clean room production	
	No release agent is used	
Classification	GMP Compliant	
Dimensions	 <p>The drawing shows a rectangular plate with rounded corners. The overall length is 127.3 units and the overall height is 85.2 units. A detail on the right shows the well depth is 21.6 units.</p>	
Colour	Clear	
Nominal Capacity	12-well	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 1	
	bag per case: 100	
Storage	Storage Temperature: 20 ° C-25 ° C, humidity 40% -60%	
	Keep clean, ventilated and out of direct sunlight	