


Technical Datasheet

Article - No:	TDS-23-022	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-PT096, 96-well, Cell Culture Plate, TC Treated	
Production conditions	Class 100,00 clean room production	
	No release agent is used	
Classification	GMP Compliant	
Dimensions		
Colour	Clear	
Nominal Capacity	96-well	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 1	
	bag per case: 100	
Storage	Storage Temperature: 20 ° C-25 ° C, humidity 40% -60%	
	Keep clean, ventilated and out of direct sunlight	