

Technical Datasheet

Article - No:	PG-QA-R025-01-0053	Revision: A0
Remarks	single use only	
Specification		
Product name	UFPT-F-200ER-L 200uL filter tips Low retention	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p>Technical drawing of a filter tip showing dimensions A through L. A is the inlet diameter, B is the outlet diameter, C is the total length, D is the length of the filter zone, E is the length of the filter zone, F is the diameter of the filter zone, G is the diameter of the filter zone, H is the diameter of the inlet, I is the diameter of the inlet, L is the diameter of the filter zone.</p>	
Colour	clear	
Graduation	moulded rings	
Nominal Capacity	I-200UL	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per case: 50	
Storage	Storage Temperature: 2 ° C-40 ° C, humidity 40% -80%	
	Keep clean, ventilated and out of direct sunlight	