


## Technical Datasheet

<b>Article - No:</b>	PG-QA-R025-01-0090	Revision: A0
<b>Remarks</b>	single use only	
<b>Specification</b>		
<b>Product name</b>	UFPT-REF-10ET-L, 10uL Extra long Pipette tips, Refull pack, Low retention	
<b>Production conditions</b>	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
<b>Classification</b>	CE/IVD	
<b>Dimensions</b>	 <p>The drawing shows a side view of a pipette tip. The total length is 45.7 ± 0.2 mm. The diameter at the narrow tip is 0.8 ± 0.02 mm. The diameter at the wider base is 6.4 ± 0.02 mm. The tip has a tapered shape with several rings near the base.</p>	
<b>Colour</b>	Clear	
<b>Graduation</b>	moulded rings	
<b>Nominal Capacity</b>	0.1-10UL	
<b>Sterile</b>	Yes	
<b>Material Characteristics</b>	polypropylene(PP),premium surface	
	pack:Polyethylene terephthalate(PET)	
<b>Product Characteristics</b>	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
<b>Shelf - Life</b>	36 months at controlled storage conditions	
<b>Packaging</b>	pcs per rack: 96	
	rack per pack: 9	
	pack per case: 10	
	pcs per case: 8640	
<b>Storage</b>	Storage Temperature: 2 ° C-40 ° C, humidity 40% -80%	
	Keep clean, ventilated and out of direct sunlight	