


Technical Datasheet

Article - No:	PG-QA-R025-01-0089	Revision: A0
Remarks	single use only	
Specification		
Product name	UFPT-REF-10ET, 10uL Extra long Pipette tips, Refull pack	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	0.1-10UL	
Sterile	Yes	
Material Characteristics	polypropylene(PP),premium surface	
	pack:Polyethylene terephthalate(PET)	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogens	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per pack: 9	
	pack per case: 10	
	pcs per case: 8640	
Storage	Storage Temperature: 2 ° C-40 ° C, humidity 40% -80%	
	Keep clean, ventilated and out of direct sunlight	